

Serravalle, 04/02/2021

SPLOŠNA BOLNIŠNICA NOVO MESTO
ŠMIHELKA CESTA 1
SI-8000 NOVO MESTO

Declaration for Uninail femoral nails and related screws – MRI Compatibility

The undersigned company LSM-MED S.r.l. based in San Marino, Strada Borrana 38 Serravalle, declares that Uninail femoral nails and the screws used for implanting these nails are manufactured in titanium alloy Ti6Al4V (ISO 5832/3).

The aforementioned material is not susceptible to magnetic resonance field; therefore, this caused neither movement, deflection, overheating nor image artifacts.

Reports

The hospitals where equipments are used for magnetic resonance are responsible for all aspects of the system which are related to the system itself: location, use, people who are subjected to MRI and that are involved in the use of the implant.

A patient who was implanted with a LSM-MED S.r.l. device composed of titanium alloy can be subjected to MRI with no risk.

However, LSM-MED S.r.l. assumes no responsibility for any damage caused by the MRI.

LSM-MED S.r.l. recommends, before use, to consult with the experts and manufacturers of the equipments for magnetic resonance imaging to learn about the compatibility of the system.

Bibliography

1. High-Field-Strength MR Imaging and Metallic Biomedical Implants: An Ex Vivo Evaluation of Deflection Forces
Frank G. Shellock, John V. Crues - AJR 151:389-392, August 1988
2. Biomedical Implants and Devices: Assessment of Magnetic Field Interactions With a 3.0-Tesla MR System
Frank G. Shellock, PhD - Journal Of Magnetic Resonance Imaging 16:721-732 (2002)
3. MR Procedures and Biomedical Implants, Materials, and Devices: 1993 Update
Frank G. Shellock, PhD, Stacy Morisoli, Emanuel Kanal, MD
4. MEDICAL DEVICES AGENCY - Guidelines for Magnetic Resonance Equipment in Clinical Use
5. MR Safety
Carolyn Roth, R.T.(R)(M)(MR)(CT)(CV)
6. Quantification of susceptibility artifacts produced on high-field magnetic resonance images by various biomaterials used for neurosurgical implants. Technical note.
Matsuura H et al. - J Neurosurg. 2002 Dec;97(6):1472-5.
7. ASTM F2119-01, Standard Test Method for Evaluation of MR Image Artifacts From Passive Implants
8. <http://www.mr-tip.com/serv1.php?type=db1&db=Artifact>
9. <http://www.radiology.upmc.edu/MRsafety/>
10. <http://www.mrisafety.com>
11. <http://www.imrser.org/>

Sincerely,

LSM-MED S.r.l.
Strada Borrana, 38
47899 Serravalle - Rep. San Marino
C.O.F. SM 27765
Iscri. Reg. Soc. 7984 del 28/06/2018

Massimiliano Manobianco
LSM-MED C.E.O.